

## EXPERIENCE MAXIMUM STABILIZATION

PRECISION SPINE  
**REFORM<sup>®</sup>MC**  
MIDLINE CORTICAL SCREW SYSTEM

Low profile 4.75 & 5.5mm cobalt chrome tulips  
**REDUCE ANATOMICAL INTERFERENCE WITHOUT COMPROMISING STRENGTH**

Audible click & consistent attachment force  
**ENSURES THAT THE TULIP IS CONNECTED**

6mm dual diameter stepped thread  
**OPTIMIZES CORTICAL BONE PURCHASE DURING FINAL SEATING**



Modular design  
**MAXIMIZES DIRECT VISUALIZATION THROUGH A LESS INVASIVE INCISION**

Enlarged T25 drive feature  
**OFFERS STABILITY DURING INSERTION**

Dual lead & cortical cancellous thread  
**PROVIDES COMPATIBLE FIXATION OPTIONS FOR PLIF & TPLIF APPLICATIONS**

Self starting aggressive distal screw tip  
**OPTIMIZES INITIAL BONE ENGAGEMENT**



**PRECISION SPINE<sup>®</sup>**  
*Discover the Difference*



## SYSTEM TRAYS

P/N	Description	P/N	Description
59-BK-1011	Reform MC 4.75 Implants (4.5-6.5mm)	39-BK-0811	Reform MC 5.5 Implants (5.5-7.5mm)
59-BK-1012	Reform MC 4.75 Implants (4.5-7.5mm)	39-BK-0812	Reform MC 5.5 Implants (4.5-9.5mm)
59-BK-2000	Reform MC 4.75 Instruments Tray 1	39-BK-0800	Reform MC 5.5 Instruments Tray 1
59-BK-3000	Reform MC 4.75 Instruments Tray 2	39-BK-0502	Reform MC 5.5 Instruments Tray 2

## BONE SCREWS

### Ø 4.5mm SCREWS

P/N	Description
59-BP-4525	Reform MC Screw - 4.5 x 25mm
59-BP-4530	Reform MC Screw - 4.5 x 30mm
59-BP-4535	Reform MC Screw - 4.5 x 35mm
59-BP-4540	Reform MC Screw - 4.5 x 40mm
59-BP-4545	Reform MC Screw - 4.5 x 45mm

### Ø 5.5mm SCREWS

P/N	Description
59-BP-5525	Reform MC Screw - 5.5 x 25mm
59-BP-5530	Reform MC Screw - 5.5 x 30mm
59-BP-5535	Reform MC Screw - 5.5 x 35mm
59-BP-5540	Reform MC Screw - 5.5 x 40mm
59-BP-5545	Reform MC Screw - 5.5 x 45mm
59-BP-5550	Reform MC Screw - 5.5 x 50mm
59-BP-5555	Reform MC Screw - 5.5 x 55mm

### Ø 7.5mm SCREWS

P/N	Description
59-BP-7530	Reform MC Screw - 7.5 x 30mm
59-BP-7535	Reform MC Screw - 7.5 x 35mm
59-BP-7540	Reform MC Screw - 7.5 x 40mm
59-BP-7545	Reform MC Screw - 7.5 x 45mm
59-BP-7550	Reform MC Screw - 7.5 x 50mm
59-BP-7555	Reform MC Screw - 7.5 x 55mm
59-BP-7560	Reform MC Screw - 7.5 x 60mm
59-BP-7570	Reform MC Screw - 7.5 x 70mm
59-BP-7580	Reform MC Screw - 7.5 x 80mm

### Ø 5.0mm SCREWS

P/N	Description
59-BP-5025	Reform MC Screw - 5.0 x 25mm
59-BP-5030	Reform MC Screw - 5.0 x 30mm
59-BP-5035	Reform MC Screw - 5.0 x 35mm
59-BP-5040	Reform MC Screw - 5.0 x 40mm
59-BP-5045	Reform MC Screw - 5.0 x 45mm

### Ø 6.5mm SCREWS

P/N	Description
59-BP-6525	Reform MC Screw - 6.5 x 25mm
59-BP-6530	Reform MC Screw - 6.5 x 30mm
59-BP-6535	Reform MC Screw - 6.5 x 35mm
59-BP-6540	Reform MC Screw - 6.5 x 40mm
59-BP-6545	Reform MC Screw - 6.5 x 45mm
59-BP-6550	Reform MC Screw - 6.5 x 50mm
59-BP-6555	Reform MC Screw - 6.5 x 55mm
59-BP-6560	Reform MC Screw - 6.5 x 60mm

### Ø 8.5mm & 9.5mm SCREWS

P/N	Description
59-BP-8540	Reform MC Screw - 8.5 x 40mm
59-BP-8545	Reform MC Screw - 8.5 x 45mm
59-BP-8550	Reform MC Screw - 8.5 x 50mm
59-BP-8555	Reform MC Screw - 8.5 x 55mm
59-BP-8560	Reform MC Screw - 8.5 x 60mm
59-BP-8570	Reform MC Screw - 8.5 x 70mm
59-BP-8580	Reform MC Screw - 8.5 x 80mm
59-BP-9560	Reform MC Screw - 9.5 x 60mm
59-BP-9570	Reform MC Screw - 9.5 x 70mm
59-BP-9580	Reform MC Screw - 9.5 x 80mm



## LOCKING CAPS, CONNECTORS, DOMINOES & RODS

### TULIPS

P/N	Description
59-MT-0301	Standard Tulip - 4.75mm
39-MT-0301	Standard Tulip - 5.5 mm
39-MT-0302	Reduction Tulip - 5.5 mm

### LOCKING CAP

P/N	Description
59-LS-0100	4.75mm Lock Screw (Gold)
39-LS-0100	5.5mm Lock Screw (Silver)



### Ø 4.75mm ROD-TO-ROD CONNECTORS

P/N	Description
59-CC-0030	Cross-Connector - 30-34mm
59-CC-0033	Cross-Connector - 33-39mm
59-CC-0038	Cross-Connector - 38-49mm
59-CC-0047	Cross-Connector - 47-67mm



### Ø 5.5-4.75mm PARALLEL DOMINOES

P/N	Description
59-DA-0203	5.5-4.75mm Open-Closed, Narrow
59-DA-0207	5.5-4.75mm Open-Closed, Wide



### Ø 4.75mm STRAIGHT CoCr RODS

P/N	Description
59-SC-5200	CoCr Straight Rod - 4.75 x 200mm
59-SC-5400	CoCr Straight Rod - 4.75 x 400mm

### Ø 4.75mm LORDOTIC CoCr RODS

P/N	Description
59-LC-5035	CoCr Lordotic Rod - 4.75 x 35mm
59-LC-5040	CoCr Lordotic Rod - 4.75 x 40mm
59-LC-5045	CoCr Lordotic Rod - 4.75 x 45mm
59-LC-5050	CoCr Lordotic Rod - 4.75 x 50mm
59-LC-5055	CoCr Lordotic Rod - 4.75 x 55mm
59-LC-5060	CoCr Lordotic Rod - 4.75 x 60mm
59-LC-5065	CoCr Lordotic Rod - 4.75 x 65mm
59-LC-5070	CoCr Lordotic Rod - 4.75 x 70mm
59-LC-5075	CoCr Lordotic Rod - 4.75 x 75mm
59-LC-5080	CoCr Lordotic Rod - 4.75 x 80mm
59-LC-5090	CoCr Lordotic Rod - 4.75 x 90mm
59-LC-5100	CoCr Lordotic Rod - 4.75 x 100mm
59-LC-5110	CoCr Lordotic Rod - 4.75 x 110mm
59-LC-5120	CoCr Lordotic Rod - 4.75 x 120mm



### Ø 5.5mm RODS AND CONNECTORS

P/N	Description
39-LT-XXXX	Ti Lordotic Rods
39-ST-XXXX	Ti Straight Rods
39-CC-XXXX	Cross Connectors

PRECISION SPINE®

www.precisionspineinc.com

