

EXPERIENCE FORWARD-THINKING MIS

120mm Integrated Extension Tab with low-profile outer diameter **ALLOWS SMALLER INCISION TO MINIMIZE TISSUE DISRUPTION**

Transitional Rod Slot design **SIMPLIFIES ROD INSERTION**

20mm of Internal Reduction Threads — with additional 30mm via instrumentation — **PROVIDES CONTROLLED ROD REDUCTION**

¾ closed Extended Tab Tulip **PROVIDES STABILITY DURING INSERTION/MANIPULATION**



True Minimally Invasive Instrumentation

- **MIS Instrumentation** connects to the Extended Tab Tulip, internally or externally, to prevent increasing size of the surgical incision
- **Compressors and Distractors** provide powerful and precise results directly at the screw-rod interface without increasing the skin incision
- **Internal Rod Reducer** allows total 50mm of controlled rod reduction
- **Rescue Towers** provide an additional level of security and confidence

PRECISION SPINE
REFORM[®]Ti^{MIS}CT
 MODULAR PERCUTANEOUS SCREW SYSTEM

Enlarged T25 Drive Feature **OFFERS IMPROVED STABILITY DURING INSERTION**

Triple lead thread **INCREASES SPEED & EFFICIENCY OF SCREW DELIVERY AND REDUCES FATIGUE**

Minimally Tapered Thread **LOWERS SCREW INSERTION TORQUE WHILE OPTIMIZING PULL OUT STRENGTH**

Self-starting, Aggressive Screw Tip **OPTIMIZES INITIAL BONE ENGAGEMENT**



PRECISION SPINE[®]

Discover the Difference



SYSTEM TRAYS

P/N	Description
64-BK-0101	Reform® Ti MIS CT Implant Tray
64-BK-0201	Reform Ti MIS CT Instrument Tray
64-BK-0210	Reform Ti MIS CT Add-On Instrument Tray



REFORM Ti CANNULATED MODULAR SCREW

5.5mm SCREWS

P/N	Description
39-SK-5530	5.5 x 30mm - Reform Ti Can Mod Screw
39-SK-5535	5.5 x 35mm - Reform Ti Can Mod Screw
39-SK-5540	5.5 x 40mm - Reform Ti Can Mod Screw
39-SK-5545	5.5 x 45mm - Reform Ti Can Mod Screw
39-SK-5550	5.5 x 50mm - Reform Ti Can Mod Screw
39-SK-5555	5.5 x 55mm - Reform Ti Can Mod Screw

7.5mm SCREWS

P/N	Description
39-SK-7535	7.5 x 35mm - Reform Ti Can Mod Screw
39-SK-7540	7.5 x 40mm - Reform Ti Can Mod Screw
39-SK-7545	7.5 x 45mm - Reform Ti Can Mod Screw
39-SK-7550	7.5 x 50mm - Reform Ti Can Mod Screw
39-SK-7555	7.5 x 55mm - Reform Ti Can Mod Screw
39-SK-7560	7.5 x 60mm - Reform Ti Can Mod Screw

6.5mm SCREWS

P/N	Description
39-SK-6530	6.5 x 30mm - Reform Ti Can Mod Screw
39-SK-6535	6.5 x 35mm - Reform Ti Can Mod Screw
39-SK-6540	6.5 x 40mm - Reform Ti Can Mod Screw
39-SK-6545	6.5 x 45mm - Reform Ti Can Mod Screw
39-SK-6550	6.5 x 50mm - Reform Ti Can Mod Screw
39-SK-6555	6.5 x 55mm - Reform Ti Can Mod Screw

8.5mm SCREWS

P/N	Description
39-SK-8535	8.5 x 35mm - Reform Ti Can Mod Screw
39-SK-8540	8.5 x 40mm - Reform Ti Can Mod Screw
39-SK-8545	8.5 x 45mm - Reform Ti Can Mod Screw
39-SK-8550	8.5 x 50mm - Reform Ti Can Mod Screw
39-SK-8555	8.5 x 55mm - Reform Ti Can Mod Screw
39-SK-8560	8.5 x 60mm - Reform Ti Can Mod Screw

NOTE: 9.5mm screws are available by request from Customer Service

PRECISION SPINE REFORM[®] Ti MIS CT MODULAR PERCUTANEOUS SCREW SYSTEM

LOCKING CAP, TULIP & RODS

LOCKING CAP

P/N	Description
64-LS-0100	Cap Screw, T30 Hexalobe



TULIP

P/N	Description
64-MT-0403	MIS Modular Ext. Tab Tulip



LORDOTIC Ti RODS

P/N	Description
63-LT-5035	Lordotic Ti Rod, x 35mm
63-LT-5040	Lordotic Ti Rod, x 40mm
63-LT-5045	Lordotic Ti Rod, x 45mm
63-LT-5050	Lordotic Ti Rod, x 50mm
63-LT-5055	Lordotic Ti Rod, x 55mm
63-LT-5060	Lordotic Ti Rod, x 60mm
63-LT-5065	Lordotic Ti Rod, x 65mm
63-LT-5070	Lordotic Ti Rod, x 70mm
63-LT-5075	Lordotic Ti Rod, x 75mm
63-LT-5080	Lordotic Ti Rod, x 80mm
63-LT-5090	Lordotic Ti Rod, x 90mm
63-LT-5100	Lordotic Ti Rod, x 100mm
63-LT-5110	Lordotic Ti Rod, x 110mm
63-LT-5120	Lordotic Ti Rod, x 120mm

NOTE: Additional rod sizes available by request from Customer Service



STRAIGHT Ti RODS

P/N	Description
63-ST-5400	Straight Ti Rod, x 400mm

NOTE: Additional rod sizes available by request from Customer Service



PRECISION SPINE[®]

www.precisionspineinc.com

